Characterization of a microchannel plate photomultiplier tube with a high sensitivity GaAs photocathode

John Martin, Paul Hink, Charlie Tomasetti, Joseph Wright

BURLE INDUSTRIES, INC.

13th Annual International Laser Ranging Workshop

13th Annual International Laser Ranging Workshop



MCP-PMT Developments


▶ 18mm GaAs MCP-PMT

- Very high, broad spectral response with excellent red sensitivity.
- Very good photon counting properties.
- Applications include Lidar, flourescence, chemiluminescence, ...

≻ 2" Square PlanaconTM

- Large area MCP-PMT with good photon counting properties.
- Ideal for tiling applications and good platform for multi-anode configurations.





85104 MCP-PMT

▶ 18 mm dual MCP.


- Very high cathode sensitivity (>1200 μA/lm).
- Excellent red response (>20% QE at 850nm).
- Excellent photon counting properties (peak-valley > 2:1).





Spectral Response

> Very high cathode response, >1200 µA/lm luminous sensitivity, 1500 µA/lm typical. ➢ Excellent QE ■ ~30% at 600nm ■ > 20% at 850nm Cathode ■ > 20% at 370nm




Typical Spectral Response

8 October 2002

Gain Curve

High gain devices,
 typically
 4x10⁶.






Single Electron Spectrum

At tube Gain > 5x10⁵
➢ Peak:Valley > 2:1
➢ Dark counts < 200 cps.


At tube Gain > 1x10⁶ ➤ Typical SER < 150% (SER= FWHM/Peak)





Single Electron Spectrum


>At higher gains ability to distinguish 1pe and 2pe peaks.



13th Annual International Laser Ranging Workshop

Dark Counts

- ≻ ~10,000 cps @ Room Temp.
- < 200 cps at -30°C for operating voltage that gives > 2:1 peak-valley and gain of > 5x10⁵.
- Currently limited by "front-end" noise, not the cathode.
- Slight increase in noise with increasing voltage.








Anode Pulse – Single PE

Anode Rise-Time: < 200ps (for single electrons) Anode Fall-Time: ~ 1ns


(Currently limited by anode design. Redesign is in progress.)





Anode Pulse – Multiple PE

Anode Rise-Time: Tek ~350 ps (using 28ps rise-time laser pulse @636nm) **Anode Fall-Time: 1.5 ns** \sim ***** **Transit Time Spread** (TTS): ~ 200 ps Ch1



8 October 2002



Lifetime of 85104

Tube Operating Conditions				Anode Output vs. Extracted	
Tube Gain	1×10^{6}	Gain		Charge	
Typical Count Rate	1x10 ⁵	Cts/s		90%	
Lifetime (@50% Anode Sensitivity Loss)			e Output		
Total Extracted	7.0	Coulombs	Anode	70%	
Charge	4.4×10^{19}	Electrons		60%	
1.2E+05 hrs	14 years	1.5x10 ¹⁴ photons (@30% QE)		50% 0.0 0.5 1.0 1.5 2. Extracted Charge (C	

8 October 2002

13th Annual International Laser Ranging Workshop




2.5

Pulse Linearity

Very good pulse linearity, ~ 400mA peak pulse current for < 5% linearity

> Average anode current is limited to ~300nA maximum.






Count Rate Linearity

Ra

- Good count rate linearity. \rightarrow Within 20% for 10⁶ cps at gain of 10^6
 - (counting pulses > 1/3 pe as defined at low count rates).
- Strip Current ~ 3mA.



Gated Operation


- Gating the photocathode of the 85104 results in > 10⁷ reduction in sensitivity for white light.
- Integrated gating circuit and external drive module in development for fast gating of PMT.

<u>Goals</u>: 10ns on/off response. On times: 100^{'s}ns – 10^{'s}μs Normally OFF/Gate ON



InGaAs Spectral Response

GaAs & InGaAs Cathodes



Developments

Improved anode response

 (Goal <200ps rise/fall times).

 Expanded photocathode offering, including InGaAs, Bi-alkali, and Multi-alkali.
 Custom anode designs will be considered.
 Gating (ongoing development!).



Burle Accessories

HV Power Supply <u>PF1056</u> New November 2002. Low Noise (~15mV ripple). Neg. 0-5kV. 667uA output current. Analog programmable. Inherent short circuit protection.





Comparison with Conventional PMT

Characteristic	MCP-PMT	Conventional
Construction	Simple Multiplier	Complex multiplier
	Complex Body	Simple Body
Processing	Slow & Complex	Fast & Simple
	"Clean"	"high alkali"
Effective QE	Moderate collection efficiency	<i>Good – Excellent Collection efficiency</i>
	> 20% (370 – 870nm)	25 – 33% peak QE



Comparison with Conventional PMT

Characteristic	MCP-PMT	Conventional
Timing Resolution	Excellent, risetime < 200ps	Poor – Very good
Open Area	Poor – Moderate with existing designs	Good - Excellent
Linearity	Limited average DC current (<1µA)	Good–Excellent for DC
	Excellent for pulsed light	Moderate – very good for pulsed light

